


Chair & Stand Unit


- ◆ An automatic and full functional Chair & Stand Unit, automatic functions are realized in different parts: automatic table moving, automatic tableboard up and down, automatic column up and down, automatic phoropter arm.
- ◆ European smart design style in appearance, which looks simple and fashionable.



Tablet Control System: besides operation panel, it can be controlled by Tablet.



Automatic moving phoropter arm: left and right.



Automatic table moving: left and right, up and down.



Specification

Unit size:	1200mm x 1300mm x 1800mm
Table size:	850mm x 450mm
Small table size:	620mm x 400mm
Drawer size:	520mm x 350mm
Control mode:	film button / phone remote control / flat panel remote control
Configuration equipment:	computer refractometer / electric cattle eye / slit lamp / focus meter

Table rotation angle:	90 °
Countertop left and right movement distance:	400mm
Countertop up and down distance:	150mm
Column (cattle eye arm) left and right moving distance:	400mm
Column (cattle eye arm) up and down distance:	120mm
Chair left and right movement distance:	150mm